


## SUMMARY:

One dimensional proton ( $^1\text{H}$ ) NMR spectra of biological mixtures (i.e. serum, plasma, CSF, urine) can easily be analyzed using Chenomx's NMR Suite. The Profiler module, coupled with a sophisticated and flexible *multi*-frequency compound library database, allows for spectral deconvolution, identification and quantification of single metabolites.

Chenomx has developed an integrated software suite and over 330 detailed signatures of common compounds used for metabolomics research. Spectral models of reference compounds can be searched and overlaid on the experimental spectra to interactively validate or invalidate potential metabolites. However, burdensome external web searches based on observed peak positions are required to search for metabolites not currently in the library.

In an ongoing effort to simplify the identification of metabolites directly in Profiler in a convivial and effortless way, Chenomx has begun a compound library expansion program using public NMR compound databases. An automated approach combining Chenomx's existing algorithms and technologies has been developed and tested on the Human Metabolome Database (www.hmdb.ca). Focusing strictly on spectra acquired in water, about 335 new compound signatures were added, effectively doubling the size of available compounds in Chenomx's library.



## EXAMPLE OF AUTOMATIC COMPOUND SIGNATURE CREATION FOR ARABINOSE

**TOP:** Chenomx native signature  
**BOTTOM:** HMDB auto-generated signature

### HMDB entries:

- High fidelity spectral reconstruction using pure lorentzian peak shapes.
- Line widths automatically adjusted to magnet shimming conditions.
- Peak clustering based on peak proximity and agnostic of chemical structure (*may be different than Chenomx*).
- Single pH (*vast majority at neutral pH*).
- Highly dependent on raw data quality (fid) and reported peaks.
- Absolute quantification not tested by Chenomx.

