


NMR Spectra Acquisition using METNOESY Pulse		
N007		
Version: 1.0	Date: 26-Sep-2018	
Version: 1.1	Date: 01-Oct-2018	

## NMR Spectra Acquisition using METNOESY Pulse Sequence

**Purpose** This procedure outlines the parameters for NMR data collection.

**Procedure** Using Chenomx NMR Suite, you can quantitatively analyze spectra produced by NOESY pulse sequences with  $t_{\text{mix}}$  of 100 ms. You can use other pulse sequences for identification purposes, but absolute quantification based on alternate pulse sequences may be unreliable.

### Required NMR Parameters

1D-1H NMR pulse sequence used for metabolomics samples. For the best results from Chenomx NMR Suite when analyzing your data, use the following parameters for NMR data collection:

- Temperature: 25 °C
- Acquisition time: 4 s
- Mixing time: 0.1 s
- Initial delay: 0.01 s
- Pre-saturation delay: 0.99 s
- Transients/scans: 32
- Steady state scans: 4
- Spectral width: 12 ppm

SOP prepared by: Nazanin Assempour  
 SOP verified by: Pascal Mercier  
 SOP approved by: Neil Taylor

